

Supervised By

Dr. Ibrahim Zaghloul

Dr. Desouky

TA. Mahmoud Haddad

Implemented By

Khaled Mohamed Soliman

Ziad Ayman Mohamed

Anas Barakat Ahmed

Mahmoud Ahmed Mohamed


Graduation Project Academic

Year 2022-2023

Final Presentation


OpticAI

Revolutionizing eye care with AI






Content

1/




01 Project Idea



Project Idea






Idea: diagnosis of eye diseases such as cataract, glaucoma, diabetic retinopathy, AMD and Myopia through fundus images using powered AI.







Domain: Bioinformatics , medical image diagnosis.

02 Problem Significance




What Is The Problem

We Aim to catch up with the eye disease diagnosis before it turns into vision loss, using high accuracy artificial intelligence models




Why

- Eye Diseases Can Cause Gradual Vision Loss Without Noticeable Symptoms.
- Most Of These Diseases Cause Irreversible Damage To The Eyes.
- Manual Examinations Are Time-Consuming, Subjective, And Can Lead To Incorrect Diagnoses And Treatments.
- The Increasing Number Of People Suffering From Eye Diseases, Particularly In Developing Countries.



How





Dataset (Fundus)


Fundus images are photographs of the back of the eye, including the retina, optic disc, macula, and fovea.

These images are captured using specialized cameras

What does our eye look like when the doctor has a look?




Models




Cataract
ResNet50
99%

1058




Glaucoma
VGG19
99%

1022



DR
InceptionV3
95%

2800




Dataset

Color Fundus Photographs/Images
From Left And Right Eyes And
Doctors' Diagnostic Keywords From
Doctors.

1051

AMD
VGG19
84%



937

Mypoia
VGG19
97%



620

Hypertension
VGG19, EfficientNetB0
72%

Cataract Assessment

91%

Our old work
CNN efficientBO
Cataract 400 fundus images

96%

Related work
CNN
Cataract ODIR-5k 1058
fundus images


99%

Final



Cataract

Dataset


Ocular Disease Intelligent Recognition (ODIR)
• 1058 Fundus Images



Preprocessing



Model Architecture



Glaucoma

Assessment

Final

83% 93%

99%

Our old work

CNN 12 layer

Glaucoma 1253 fundus
images

Related work

CNN


LAG database
fundus images cropped to
optic cup

Glaucoma

Dataset


GlaucomaDataset From Kaggle

- 1022 Fundus Images




Confusion Matrix


99% acc



Preprocessing



Model Architecture



Diabetic Retinopathy

Assessment

Final

75%

Our old work
CNN 7 layers
Messidor-1 1744 images

96%

Related work
unknown model
Messidor-2


95%

Diabetic Retinopathy


Dataset

Asia Pacific Tele-Ophthalmology Society(Aptos-2019-Blindness)


- 2800 Fundus Images




Preprocessing



Confusion Matrix **95% acc**




Model Architecture



Experts Feedback

(Dr. Amro Al Shafi – PHD In Ophthalmology)



1- Common Diseases

- age-related macular degeneration (AMD)
- Hypertension
- Myopia

2- Simultaneously Diagnosis

we provide 5 eye classification models that can work simultaneously to diagnose multiple diseases at once.

3- Motivation

Doctor inspired us to continue on such project as it can be used for insurance and periodic check up

AMD Assessment

Final

82%

Our old work

InceptionV3

ODIR5K

91%

Related work

Nine layer CNN model

Private dataset


84%

AMD


Dataset

Ocular Disease Intelligent Recognition (ODIR5K)


- 1051 Fundus Images




Confusion Matrix **84% acc**



Preprocessing



Model Architecture



Myopia Assessment

94%

Our old work

ResNet50, efficientB7

ODIR5K

98%

Related work

VGG19

ODIR5K

97%


Final

Myopia


Dataset

Ocular Disease Intelligent Recognition (ODIR5K)


- 973 Fundus Images




Preprocessing



Confusion Matrix **97% acc**







Model Architecture







OpticAI

Mobile App







OpticAI

Mobile App



OpticAI

Mobile App



OpticAI

Feedback

The screenshot displays a web browser window for 'opticAI' at 'opticAI.com/Feedback'. The main content area shows two feedback posts. The first post, from 'example@gmail.com' (posted 'About a minute ago'), expresses satisfaction with the ML model's accuracy in identifying eye diseases. The second post, also from 'example@gmail.com' (posted 'About a minute ago'), provides a more detailed compliment, stating the model correctly identified the user's condition and recommended appropriate next steps, revolutionizing eye healthcare.

Below the posts, device information is listed:


- Locale: en-US (USA)
- Platform: Android
- SDK Version: 1.5.0 (outdated)
- App Version: not set
- Build Number: not set
- Dart Version: 2.18.6

At the bottom, there are links for 'User Email' (example@gmail.com) and a 'Copy' button.

On the right side of the browser window, a mobile application interface for 'Eye Diseases Service' is overlaid. It features a green header with a back arrow and the service name. Below it is a section for 'Upload Image' with 'From Gallery' and 'From Camera' buttons. A thumbnail image of an eye fundus photograph is shown. The results section displays the prediction: '77' and the result: 'Cataract'. A large green button labeled 'Consult Doctor' is present, along with a note: 'You could chat with doctor to consult them'.

The entire browser window is highlighted with a thick red border.

Website



Our services

Eye Diagnosis Services



99% Accuracy

99% Accuracy

95% Accuracy

Website

The screenshot shows a web browser window for the website opticAI.com. The page is titled "Our services" under "Eye Diagnosis Services". It displays five service cards:

- Cataract**: Features two images of children holding soccer balls, one labeled "Normal vision" and one labeled "Cataract". Accuracy: 99%.
- Glaucoma**: Features two images of children holding soccer balls, one labeled "Normal vision" and one labeled "Glaucoma". Accuracy: 99%.
- Diabetic Retinopathy**: Features two images of children holding soccer balls, one labeled "Normal vision" and one labeled "Diabetic retinopathy". Accuracy: 95%.
- Macular Degeneration**: Features two images of children holding soccer balls, one labeled "Normal vision" and one labeled "Macular degeneration". Accuracy: 84%.
- Myopia**: Features two images of children holding soccer balls, one labeled "Normal vision" and one labeled "Myopia". Accuracy: 96.5%.

A yellow callout bubble with the text "Click Here" points to the "Start Service" button of the Cataract card. A yellow arrow points from the "Click Here" text towards the "Start Service" button.

Our services

Eye Diagnosis Services

 
Normal vision Cataract

 
Normal vision Glaucoma

 
Normal vision Diabetic retinopathy

 
Normal vision Macular degeneration

 
Normal vision Myopia

Cataract

DenseNet-201 is a deep learning model architecture

eye condition characterized by the clouding of the natural lens in the eye

Start Service

Glaucoma

VGG19 consists of 19 layers pretrained model in computer vision tasks

Glaucoma is a leading cause of irreversible blindness worldwide

Start Service

Diabetic Retinopathy

Inception-v3 architecture consists of 48 layers

eye condition that affects people with diabetes

Start Service

Macular Degeneration

Myopia

Click Here

Website

The screenshot shows a web browser window for 'opticAI' at 'opticAI.com'. The main heading is 'Cataract Service Model'. A subtext states: 'Our Cataract Service Model is designed to provide comprehensive care and support for individuals diagnosed with cataracts, ensuring a streamlined and personalized experience from diagnosis to treatment.' Below this are two circular progress indicators: one green circle labeled '99%' for 'Model Accuracy' and one grey circle labeled '2%' for 'Real-life error'. A call-to-action button says 'Start service by Uploading fundus Image' with a placeholder 'Click to upload PNG, JPG (MAX. 250 X 250 px)'. Below this is a video thumbnail titled 'Animation: Cataract' showing a cross-section of an eye with a label 'Retina' pointing to the inner wall of the eye.

Cataract Service Model

Our Cataract Service Model is designed to provide comprehensive care and support for individuals diagnosed with cataracts, ensuring a streamlined and personalized experience from diagnosis to treatment.

Model Accuracy 99%

Real-life error 2%

Start service by
Uploading fundus Image


Click to upload
PNG, JPG (MAX. 250 X 250 px)

Learn more about Cataract

Animation: Cataract

Retina

OpticAI API




```
● ● ●

1 from flask import Flask, request, jsonify
2 from tensorflow import keras
3
4 from keras_preprocessing.image import load_img
5 from keras_preprocessing.image import img_to_array
6 import numpy as np
7 from PIL import Image
8 from flask_cors import CORS
9 import cv2
10 import numpy as np
11
12
13 app = Flask(__name__)
14 CORS(app)
15
16 cataract_classification_model = keras.models.load_model(
17     "D:/Ocular/Cattract/OcularCattract/my_model/ResNet50-Cataract_classification-98.98.h5")
18 glaucoma_classification_model = keras.models.load_model(
19     "D:/Ocular/Cattract/OcularCattract/my_model/VGG19-Glaucoma_classification-98.91.h5")
20 retinopathy_classification_model = keras.models.load_model(
21     "D:/Ocular/Cattract/OcularCattract/my_model/inception-DR_classification-95.00.h5")
22 amd_classification_model = keras.models.load_model(
23     "D:/Ocular/Cattract/OcularCattract/my_model/VGG19-AMD_classification-84.21.h5")
24 myopia_classification_model = keras.models.load_model(
25     "D:/Ocular/Cattract/OcularCattract/my_model/VGG19-Myopia_classification-96.59.h5")
```

OpticAI

API



Postman screenshot showing a POST request to http://localhost:4000/cataract with a file parameter named 'imagefile' containing 'cataract_011.png'. The response body is a JSON object with keys: disease, prediction, and result.

Request Details:

- Method: POST
- URL: http://localhost:4000/cataract
- Body Type: form-data
- Body Parameters:

Key	Value	Description
imagefile	cataract_011.png	
Key	Value	Description


Response Headers:

- Status: 200 OK
- Time: 2.16 s
- Size: 255 B

Response Body (Pretty JSON):

```
1  "disease": "Cataract",
2  "prediction": 87,
3  "result": "Negative"
```

Pros & Cons





Positives


- ✓ cover a set of models that diagnose various eye diseases with a higher level of accuracy
- ✓ communication with various ophthalmologists

Negatives

- ✗ Requires special fundus cameras which is relatively expensive (3000\$ -> 7000\$)
- ✗ Input Data might not be enough for certain diseases



03 Project Specifications



Functional & Non-Functional

Functional


- Sign in / Sign Up
- Image Upload
- Machine Learning Model
- User Feedback
- Chat Functionality

Non-functional


- Security
- Usability
- Reliability
- Compatibility

Stakeholders


- Patients
- Ophthalmologists
- Researchers
- Medical Equipment Manufacturers




System Architecture




Class Diagram




Use Case Diagram






Entity Relationship Diagram






Sequence Diagram






04 Time Plan








Time Plan



05 Conclusion






Conclusion

- The proposed solution will be integrated into a user-friendly interface for ease of use.
- This project represents a step towards better healthcare outcomes and opportunities for further research and development in the area.



Looking Into The Future



Thank You!